



Emergency Preparedness Systems



Pediatric TEMPS Bed

The standard adult model of the original **TEMPS Bed** has been downsized for the pediatric market and includes all the standard features of the original **TEMPS Bed**. Provides immediate additional bed capacity for E.D. surge overflow or MCI. Assembles in 30 seconds and ready for patient use. The **Pediatric TEMPS Bed** is constructed of the highest quality poly-resin structural tubing which will not chip, dent, peel, fade or sustain bacterial growth. The patient is supported by a vinyl coated polyester mesh that is tension adjustable, removable and easily decontaminated.

Pediatric TEMPS Bed

Model	Description
PTB500 System	5 each PTB100, pediatric <i>TEMPS Beds</i> with triage I.D. System and 1 each PTB Cart, transport and storage cart with waterproof dust cover
PTB100	1 each Pediatric <i>TEMPS Bed</i> with triage I.D. system
PTB Cart	1 each transport and storage cart with dust cover holds 5 each PTB100's

OPTIONS:

Same options available as the standard adult model.

Pediatric TEMPS Bed Specifications:

- Overall Length: 52"
- Overall Width: 25"
- Folded dimensions: 25" x 52" x 11"
- Working height: 29"
- Weight capacity: 350 lbs