

Shelter bed

Ideal for isolation, disaster preparedness, post acute surge overflow or evacuation shelter. Includes 4 position elevating back rest, 2" poly mattress with water proof cover and 650 lb capacity. System utilizes fold, stack and store technology with a mobile transport and storage cart that holds 6 beds per system.



The Shelter Bed has been developed in conjunction with leading experts in the field of disaster preparedness and has been tested in rigorous patient overflow scenerios.



Folds flat for transport and storage and assembles in 20 seconds.



Snap on I.V. Pole




Mobile transport and storage cart with dust cover holds 6 beds.

SPECIFICATIONS:
 Mattress Length: 80"
 Mattress Width: 30"
 Weight Capacity: 650 lbs
 Bed Height
 Folded for Storage: 8"
 In Use: Customer Specific
 Overall Length: 84" in use
 Overall Length: 88" folded
 Overall Width: 36.5"
 10 year warranty.



The Shelter Bed is constructed of the highest quality polyvinyl chloride resin structural tubing that will not chip, dent, peel, stain, fade or sustain bacterial growth.



OPTIONS:
 Side Rails
 I.V. Pole
 3" Mattress
 Disposable Pillow
 O2 Holder
 Patient Chart Holder
 4" Mattress
 Disposable Sheets