

TEMPS transport stretcher system

An economical, multi-function stretcher system for MCI's, decon, equipment staging, rescue, transfer and rehab.

The TEMPS transport stretcher system is ideal for decon, rapid intervention, patient transport, equipment transfer and rehab. Constructed from the highest grade structural poly-resin, the TEMPS stretcher was developed in conjunction with leading experts and is a proven system for MCI's.



Overall width of 25" will pass through most doors.



The TEMPS stretcher doubles as an equipment transfer system and shortens staging time and manpower requirements.

TEMPS Stretcher specs:

Overall Length: 75"

Overall Width: 25"

Weight capacity: 300 lbs,

Extra weight capacity models available.

Patent Pending



TEMPS Bed mesh sling allows water to drain through to collection point and is easily removed for cleaning.



Low center of gravity makes patient transfers an easy, one man function.



Multi-purpose stretcher is ideal for rehab.