

HEALTHMARK

Pharmacy Accessory Label Printer (PALP)

BARCODE LABELING SOLUTION FOR ANCILLARY PACKAGING

If your goal is attaining Bedside Medication Verification (BMV) then PALP accomplishes this through a Pharmacy software labeling program with specialty labels to enable barcoding and strong adhesion for the many uniquely shaped ancillary medications, such as:

- Oral syringes
- Ampoules
- Small vials
- Nebules
- Inhalers
- IV bags
- Creams
- Ophthalmics



Thermal Transfer
Barcode Printer

PALP is a turnkey solution for labeling all of your ancillary medications!

Includes Zebra GK420 printer, custom Pharmacy Pak-EDGE™ Barcode Labeling Software, specialty format and adherence labels.

TECHNICAL SPECIFICATIONS

Dimensions	8.25"L x 6.75"W x 6.0"H
Weight	3 lbs.
Power	120V, 50/60Hz

Pak-EDGE™ Custom Pharmacy UD Barcode Label Software

- User-maintained product master file for quick recall of standard labels
- Enhanced security options allow control of data files and labeling options
- **Supports all Linear, 2D, and GS1 Barcode Formats**
- Live label preview
- Pharmacy Software – Lot Recall, Expiry Dating, Tall Man Lettering, Reports and much more.
- Templates accommodate label printing in conjunction with ISMP guidelines



Pink Circle in Square Label
(GA-CDP) on **Patch**




Clear Flag Label
(GA-FL) on
Neomed® Amber Oral Syringe




Paper Rectangular Label
(GA-L5) on **IV Syringe**



Side Perf Label(GA-SP) on
Neomed® Clear Oral Syringe



Cryogenic Rectangular Label (GA-L2) on Vial



Paper Rectangular Label (GA-L5) on TampAlert®



Circle Label (GA-C) on Inhaler



Cryogenic Rectangular Label (GA-L1) on Nebule



Paper Rectangular Label (GA-L5) on 250 mL IV Bag



Circle in Square (GA-CD) on Ampoule

Specialty Labels to meet product needs

Side Perf: Easy application and strong adhesion for oral syringes

Clear Flag: Labels with clear stem – print on one or both sides

Cryogenic:


- Smudge proof with extra adhesive
- For fridge/freezer applications

Rectangular: For small to large items including IV bags

Circle: Small size bar code only application with ability for minimal written text such as: bar code, lot #, and/or expiry date

Circle in Square:

- For small size bar code only application
- Two colour options available (white and pink)



Side Perf




Clear Flag Label



Circle in Square



Circle



Cryogenic



Rectangular

Ask your local Healthmark Rep for complete listing of label options.

1.800.665.5492 · info@healthmark.ca

www.healthmark.ca