

Finally, an automated solution for filling oral liquid syringes!



Introducing the Precifill™ Dispensing Pump

Another medication safety solution from Healthmark!

- Simple to use
- Significantly faster than manually filling oral liquid syringes
- Reduces margin of error – accurately dispenses medications
- Improves time management for pharmacy staff



- Reduces repetitive strain of manual filling
- Small footprint supports organized workflow
- User-friendly software
- Flexible – can also be used for IV preparations

HEALTHMARK

Precifill

Introducing the Precifill™ Pump

Versatility and Flexibility

In addition to automated filling of oral syringes, the Precifill™ Pump is ideal for filling IV products that require a high level of accuracy, such as syringes, mini-bags, cassettes and more.

Medication Safety

Fills oral liquid syringes or IVs with high-precision accuracy $\pm 1\%$ (volumes $>1\text{ml}$) and $\pm 2\%$ (volumes $<1\text{ml}$).

Accurate filling of volumes as low as 0.2ml and syringe sizes as small as 1ml .

Eliminates manual filling and reduces the possibility of repetitive strain injuries.

Improved Pharmacy Function

The Precifill™ Pump has a unique capability to accurately fill small syringes with volumes such as oral morphine or IV oxytocin. The machine also delivers consistent and repeatable filling with viscous or non-viscous solutions.

Users will appreciate the fast filling speed – a few seconds per syringe depending on viscosity and fill volume. Best of all, it is simple to set up and use.

User-friendly Software

The intuitive touch control screen allows users to set up any production run in minutes. The one step “Adjust” calibration function and memory program storage capability provides ease-of-use in any situation. For user convenience, the software is also multilingual.



CODE	DESCRIPTION	PACK
BX-PF100	Precifill™ Pump (automated filling of oral liquid syringes and IV drugs)	
BX-3051	High Flow Set	25/cs
BX-3051-2	Low Flow Set	25/cs
BX-3051-1	Sterile Liquids Set with vented spike	25/cs

HEALTHMARK

Toll Free: (800) 665-5492
Fax: (514) 336-7111
e-mail: info@healthmark.ca
www.healthmark.ca



Scan here for
Healthmark
website