

Low Server Specifications - RXMDAP-LS1

| ITEM | DESCRIPTION |
|---|---|
| Processor (minimum) | Intel® Xeon® Bronze 3106 1.7G, 8C |
| Operating System | Microsoft™ Windows® Server 2016 Standard Edition |
| Database | Microsoft™ SQL Server® 2017 Standard OEM w/5 CALs |
| Memory (minimum) | 16 GB Memory (2) 8GB RDIMM, 2666MT/s, Single Rank |
| Storage (minimum) | RAID 5 (3) 600GB 10K RPM SAS 12Gbps 512n 2.5in Hot-plug Hard Drive |
| Network Adapter | Intel Ethernet I350 DP 1Gb Server Adapter |
| Media Drive | DVD+/-RW, SATA |
| Integrated Remote Access Controller (recommended) | 1DRAC7 Enterprise |
| Raid (recommended) | RAID5 512MB NV cache |
| Power Supply | Dual Hot-Plug Power Supply 550W |
| Form Factor | Rack Mount (1U) |
| Additional | Sliding Rails w/Cable Management Arm |
| | Power Distribution Board for Power Supplies |
| | Power cord w/wall plug (varies by region) |
| Service | Mission Critical Package (4-hour 7x24 on-site service w/emergency dispatch) |
| | 3 Years Pro Support 7x24 HW/SW Technical Support & Assistance |
| | Hardware Limited Warranty Plus Next Business Day On-Site Service |

Mid Server Specifications - RXMDAP-MS1

| ITEM | DESCRIPTION |
|---|---|
| Processor (minimum) | 2 x Intel® Xeon® Bronze 3106 1.7G, 8C |
| Operating System | Microsoft™ Windows® Server 2016 Standard Edition |
| Database | Microsoft™ SQL Server® 2017 Standard OEM w/5 CALs |
| Memory (minimum) | 64 GB Memory (8) 8GB RDIMM, 2666MT/s, Single Rank |
| Storage (minimum) | RAID 10 (4) 600GB 10K RPM SAS 12Gbps 512n 2.5in Hot-plug Hard Drive |
| Network Adapter | Intel Ethernet I350 DP 1Gb Server Adapter |
| Media Drive | DVD+/-RW, SATA |
| Integrated Remote Access Controller (recommended) | 1DRAC7 Enterprise |
| Raid (recommended) | RAID1+0 1GB NV Cache |
| Power Supply | Dual Hot-Plug Power Supply 550W |
| Form Factor | Rack Mount (1U) |
| Additional | Sliding Rails w/Cable Management Arm |
| | Power Distribution Board for Power Supplies |
| | Power cord w/wall plug (varies by region) |
| Service | Mission Critical Package (4-hour 7x24 on-site service w/emergency dispatch) |
| | 3 Years Pro Support 7x24 HW/SW Technical Support & Assistance |
| | Hardware Limited Warranty Plus Next Business Day On-Site Service |

High Server Specifications - RXMDAP-HS1

| ITEM | DESCRIPTION |
|---|---|
| Processor (minimum) | 2 x Intel® Xeon® Bronze 3106 1.7G, 8C |
| Operating System | Microsoft™ Windows® Server 2016 Standard Edition |
| Database | Microsoft™ SQL Server® 2017 Standard 4CORE OEM |
| Memory (minimum) | 64 GB Memory (8) 8GB RDIMM, 2666MT/s, Single Rank |
| Storage (minimum) | RAID 10 (4) 600GB 10K RPM SAS 12Gbps 512n 2.5in Hot-plug Hard Drive |
| Network Adapter | Qlogic™ FastLinQ 41112 Dual Port 10Gb SFP+ Adapter, Full Height |
| Media Drive | DVD+/-RW, SATA |
| Integrated Remote Access Controller (recommended) | 1DRAC7 Enterprise |
| Raid (recommended) | RAID1+0 1GB NV Cache |
| Power Supply | Dual Hot-Plug Power Supply 550W |
| Form Factor | Rack Mount (1U) |
| Additional | Sliding Rails w/Cable Management Arm |
| | Power Distribution Board for Power Supplies |
| | Power cord w/wall plug (varies by region) |
| Service | Mission Critical Package (4-hour 7x24 on-site service w/emergency dispatch) |
| | 3 Years Pro Support 7x24 HW/SW Technical Support & Assistance |
| | Hardware Limited Warranty Plus Next Business Day On-Site Service |

AccessRx MD Virtual Server - Software Kit - RXMDAP-VS

Customer is required to provision the virtual server to specifications provided by TouchPoint Medical, including operating system, database server(s), web server(s) and application server(s). Customer is responsible for all licensing associated with this provisioning. Customer is also required to provide remote access to qualified TouchPoint Medical technical support personnel for installing and/or configuring application software to run on the server. TouchPoint Medical technical support personnel will conform to reasonable security requirements of the facility and the company provide signed HIPAA compliance acknowledgment and/or limited business partnership agreements on request. TouchPoint Medical prefers and provides an integrated Axeda remote management solution at no additional charge with all Access platform products, or the Customer may provide their own solution at their own cost.

Americas & Global Headquarters

Oldsmar, FL, USA
Tel: (1) 800 947 3901

Europe Headquarters

Melsele, Belgium
Tel: (32) (3) 568 33 00

UK Headquarters

Rugby, United Kingdom
Tel: (44) (0) 844 576 1247

Americas & Global Headquarters

Oldsmar, FL, USA
Tel: (1) 800 947 3901

Australia / New Zealand / Singapore / Philippines / Thailand Headquarters

Victoria, Australia
Tel: (613) 9330 5050

China / Hong Kong / Macau, Taiwan Headquarters

Shanghai, China
Tel: 86 1361 1680 312