




**NXT - SERIES 40 – EXTRACTION HOOD – DUAL CHEMICAL MODEL**



Model	Dimension A	Work Zone (W x D x H)	Electrical
<input type="checkbox"/> NXT-340	35.375" (900)	31.5" x 21.5" x 22" (800 x 546 x 559)	3A, 115V, 50/60Hz
<input type="checkbox"/> NXT-440	47.25" (1200)	43.5" x 21.5" x 22" (1105 x 546 x 559)	3A, 115V, 50/60Hz

**STANDARD FEATURES:**

- ✓ Cold-rolled steel construction for exterior shell with durable, anti-bacterial powder coat finish.
- ✓ Stainless steel interior lining, baffle, and dished work surface
- ✓ Microprocessor with LCD display to control lights, fans.
- ✓ LED task lighting & UV germicidal light.
- ✓ Motorized vertical sliding face shield c/w interlocked with fans and UV light when in the closed position.
- ✓ UV resistant safety glass sash and side windows.
- ✓ Dual carbon filters pre-charged with media specific to hood use.
- ✓ Manual adjustable base stand for 29"-36" high work surface height range c/w casters for ease of mobility.
- ✓ Low flow alarm to monitor exhaust airflow.
- ✓ Quick change disposable pre-filters, MERV 8.
- ✓ Two cord pass-through grommets in side panels.
- ✓ 15' power cord with grounded plug.
- ✓ 1 year warranty on parts and labour

**OPTIONAL FEATURES:**

**Accessories**

- Stainless steel IV bar, complete with (5) hangers
- Stainless steel side panels in lieu of glass
- Transition adaptor canopy, 8" dia collar

1<sup>ST</sup> Carbon Filter Code:  MM-1000     MM-1355     MM-3000     MM-9000     Unknown, please contact manufacturer  
 2<sup>ND</sup> Carbon Filter Code:  MM-1000     MM-1355     MM-3000     MM-9000     Unknown, please contact manufacturer

Notes / Comments:

Customer Signature \_\_\_\_\_  
Date \_\_\_\_\_